

Long live the doubly charmed particle

May 25 2018



The LHCb detector seen in 2018 in its underground cavern. The excellent precision of this detector allowed LHCb physicists to perform detailed measurements on the doubly charmed particle they discovered only last year. Credit: M. Brice, J. Ordan/CERN

Finding a new particle is always a nice surprise, but measuring its characteristics is another story and just as important. Less than a year



after announcing the discovery of the particle going by the snappy name of Ξ_{cc}^{++} (Xi_{cc}⁺⁺), this week the LHCb collaboration announced the first measurement of its lifetime. The announcement was made during the <u>CHARM 2018 international workshop</u> in Novosibirsk in Russia: a charming moment for this doubly charmed particle.

The Ξ_{cc}^{++} particle is composed of two charm quarks and one up quark, hence it is a member of the baryon family (particles composed of three quarks). The existence of the particle was predicted by the Standard Model, the theory which describes elementary particles and the forces that bind them together. LHCb's observation came last year after several years of research. Its mass was measured to be around 3621 MeV, almost four times that of the proton (the best-known baryon), thanks to its two <u>charm quarks</u>.

The Ξ_{cc}^{++} particle is fleeting: it decays quickly into lighter particles. In fact it was through its decay into a Ac+ baryon and three lighter mesons, K-, π + and π +, that it was discovered. Since then, LHCb physicists have been carrying on an analysis to determine its lifetime with a high level of precision. The value obtained is 0.256 picoseconds (0.000000000000256 seconds), with a small degree of uncertainty. Though very small in everyday life, such an amount of time is relatively large in the realm of subatomic particles. The measured value is within the range predicted by theoretical physicists on the basis of the Standard Model, namely between 0.20 and 1.05 picoseconds.

To achieve this precise result, LHCb physicists compared the measurement of the lifetime of the Ξ_{cc}^{++} with that of another particle whose lifetime is well-known. They based their measurements on the same sample of events that led to the discovery.

Measuring the lifetime of a particle is an important step in determining its characteristics. Thanks to the abundance of heavy quarks produced by



the Large Hadron Collider (LHC) and the excellent precision of the LHCb detector, physicists will now continue their detailed measurements of the properties of this charming particle. With these types of measurements, they are gaining a better understanding of the interactions that govern the behaviour of <u>particles</u> containing <u>heavy quarks</u>.

Provided by CERN

Citation: Long live the doubly charmed particle (2018, May 25) retrieved 2 May 2024 from <u>https://phys.org/news/2018-05-doubly-charmed-particle.html</u>

This document is subject to copyright. Apart from any fair dealing for the purpose of private study or research, no part may be reproduced without the written permission. The content is provided for information purposes only.