


Marine biologists clarify how specialized cells in squid skin are able to control the animal's coloration

1 July 2015




This iridescent skin patch of the California market squid shows its dynamic iridocytes. The dark spots in the background are the pigment filled chromatophores. Credit: UCSB

Perhaps not the brightest of cephalopods, the California market squid (*Doryteuthis opalescens*) has amazing light-manipulating abilities. While this species shares the gift of camouflage with most other cuttlefish, octopus and squid in the cephalopod family, it is also capable of activating, shuttering and directing its own iridescence in multiple ways.

New research conducted in the UC Santa Barbara lab of Daniel Morse, professor emeritus in the Department of Molecular, Cellular and Developmental Biology (MCDB), demonstrates that the squid's ability to "tune" its colors is correlated with the presence of specific sequences of reflectins, proteins unique to the light-sensing tissue of [cephalopods](#). The findings appear in the *Journal of Biological Chemistry*.

While the diverse optical behaviors of cephalopod tissues have been attributed to reflectins for almost a decade, this work shows for the first time how reflectin protein subtype structure, localization, distribution and relative abundance correlate with the squid's optical output. Building on the lab's previous work, the new research details the mechanisms of the animal's tunable (adaptive) and nontunable (static) iridocytes. Specialized cells in squid skin, iridocytes—also known as iridophores—produce color via Bragg reflection, by which light is reflected in a very regular and predictable manner.



A close-up view of the iridocytes, showing their potential to tune to a broad range of colors. Credit: UCSB

"The discovery reported in this paper reveals the subtlety and power of the reflectin proteins to fine-tune the colors of living cells with a beauty that reminds us of paintings by Monet," said Morse, who is a co-author on the paper.

In previous papers, the researchers demonstrated


that the neurotransmitter acetylcholine triggers reflectin proteins to condense and expand deep pleats in the cell membrane through the chemical switch of phosphorylation—a process that turns [protein](#) enzymes on and off, altering their function and activity. This in turn tunes the brightness and color of the reflected light so that these layers or lamellae operate as a tunable Bragg reflector.

The newly published research characterizes the mechanistic action of iridocytes by characterizing their reflectin proteins. "It's a very complicated system," said lead author Daniel DeMartini, who conducted this research as a graduate student and is now a postdoc in another lab at UCSB. The results suggest that distinctive reflectin subtypes have different roles.

noted that because static iridocytes look thicker than the dynamic ones, their potential for increased membrane interface may require larger amounts of reflectin C. Location—responsive iridocytes are located on the squid's back with more passive ones on its underside—may also be a factor, one that requires further exploration to fully understand.

"Our new work into the cellular mechanisms of squid iridescence allows us to refine experiment design for reflectin-based adaptive photonic structures," DeMartini said. "In the future, this could be very useful for creating synthetic optical materials that—like squid skin—can be tuned."

Provided by University of California - Santa Barbara



The adaptive iridocytes in the skin of the California market squid are able to tune color through most of the spectrum. Credit: UCSB

"We found three major types of reflectins—A1 and A2, B and C—which in various combinations determine slightly different iridocyte function," he explained. "A-type reflectins are found in static iridocytes; B are enriched in fully tunable ones, which contain A as well. We also found reflectin C, a new type that is probably important for anchoring proteins to membranes."

Reflectin C is found in both static and tunable iridocytes; however, their amount varies. DeMartini

APA citation: Marine biologists clarify how specialized cells in squid skin are able to control the animal's coloration (2015, July 1) retrieved 29 September 2020 from <https://phys.org/news/2015-07-marine-biologists-specialized-cells-squid.html>

This document is subject to copyright. Apart from any fair dealing for the purpose of private study or research, no part may be reproduced without the written permission. The content is provided for information purposes only.